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CMA System

Instruction for use EN

- Description**
Rotary dental instruments of type 2 compliant with EN ISO 3630 standard. The materials used in the working part: Nickel-Titanium and stainless steel (for PreSequencer).
Intended for use in connection with an endodontic contra-angle that complies with the current EN ISO 1797 standard (type 1).
- Indication**
CMA system is indicated for root canal preparation during non-surgical endodontic treatment and retreatment.
- User**
Intended for use by dental professionals only.
- Contraindication**
There is no contraindication to the use of CMA products after the dental professional has decided to perform a non-surgical root canal treatment.
- Undesirable effects/ Complications**
In cases of complex canal anatomy, per-operative risks (breakage, compression, stripping, perforation, blockage) could occur and lead to a risk of infectious processes.
- Warnings/Precautions**

- The decision to use an endodontics instrument must be based on the clinical case expertise, particularly where the canal anatomy is considered to be too complex.
- Patients identified as having a risk of infectious endocarditis.
- Contains nickel and titanium: should not be used in individuals with a known allergy to these metals.
- Respect the good dental practice, in particular the use of a dental dam and gloves.
- Use in continuous rotation at the recommended speed.
- In cases of complex anatomy, the maximum torque may differ from the value recommended by the manufacturer.
- Use according to the recommended procedure (§11).
- Instruments supplied non-sterile:** Perform the cleaning and sterilization procedure **before each use**, in accordance with the instructions for use in §10. Adhere to the maximum number of uses indicated in §8.
- Keep the UDI information which only appears on the labelling until the last use.
- Do not use if doubt on the respect of storage conditions.
- In case of doubt concerning the product identification, do not use.
- Check the condition of the instrument blade before use between each canal. If the instrument is damaged or shows signs of wear, do not use it.
- Inform the manufacturer and the national regulatory authority of any serious incident relating to the instrument.

7. Clinical claims

CMA in normal use conditions.

- Clinical performance:**
- Shaping ability: adapted debris removal and cutting efficiency
 - Respect of canal anatomy
 - Geometry adapted to irrigation*
 - Maintains the apical foramen in its original position**

- Safety:**
- Limited risk of breakage
 - Limited peroperative risks*
 - Limited risk of screwing effect
- *Not applicable for the Coronal and Pre-Sequencer instruments.
**Not applicable for the Coronal instrument

8. Characteristics

Maximum recommended torque (not mandatory): see table 1 «CMA SYSTEM Characteristics».

- Number of reuses:**
- Maximum recommended 10 canals if the instrument is not visually damaged.
 - Maximum recommended 3 canals for the PreSequencer file.

	Ø	%	L	Stop	Ring	RPM	Recommended Torque (N.cm)
Coronal	25	10	17	-	○	300 - 400	2 - 3
Median	25	6	21-25	○	●	300 - 400	2 - 3
Apical1	20	4	21-25-29	○	●	300 - 400	2 - 3
Apical2	20	6	21-25-29	○	●	300 - 400	2 - 3
Apical3	30	6	25	○	●	300 - 400	2 - 3
Apical4	40	7	25	○	●	300 - 400	2 - 3
Pre sequencer	17	2	21-25	○	●	300 - 400	2 - 3

Table 1: CMA SYSTEM characteristics

9. Protocole


Figure 1: Protocol for use of the CMA SYSTEM

10. Instructions for use

- General recommendations**
- For all metal devices, the use of anticorrosion disinfecting and cleaning agents is recommended.
 - For your own safety, please wear personal protective equipment (gloves, glasses and mask).
 - Do not use cleaning or disinfecting agents containing phenol, aldehyde and alkaline composition.
 - Use cleaning or disinfecting agents compatible with anodized aluminum.
 - Always respect the instructions for use provided by the manufacturer of the products.
 - The user must ensure that the processing method used, including re-sources, materials and personnel, is appropriate and meets the applicable requirements.
 - The state of the art and national laws require that validated procedures be followed.

Material needed	<ul style="list-style-type: none"> Gloves, masks, gown as recommended by the manufacturer of cleaning agent and detergent Tap or deionized water Disinfectant (neodisher® Septo Active) Detergent (neodisher® MediZym) Small soft brushes Container Ultrasonic tub or washer-disinfector <p><i>Note: All material used should be cleaned and replaced regularly. Identify material used for each step of the process (initial treatment, cleaning or rinsing).</i></p>
------------------------	---

1	<p>Initial treatment</p> <p>Immerse the products in a container or on a wipe with tap water at 20-40°C and 1.0% neodisher® Septo Active for 5 to 15 minutes. Rinse the products with tap water at 20-40°C for 1 minute. If there is waiting time before the next step, make sure the device stays moist by placing it in a wet wipe. Do not exceed 1 hour waiting time.</p> <p><i>Note:</i></p> <ul style="list-style-type: none"> Do not use fixing agents or hot water (> 40 °C), as this causes fixation of residues and can impair successful cleaning. Follow instructions and respect concentrations and immersion times specified by the manufacturer (an excessive concentration may cause corrosion or other defects on devices).
----------	--

2	<p>Preparation before cleaning</p> <p>If the devices have visible impurities, it is recommended to manually pre-clean by brushing under tap water at 20-40°C for at least 1 min with a soft brush until all impurities have been removed.</p> <p><i>Note: Follow instructions and respect concentrations and immersion times specified by the manufacturer (an excessive concentration may cause corrosion or other defects on devices).</i></p>
----------	---

3	<p>Visual inspection</p> <p>Inspect the used products and discard damaged products (broken, untwisted or abnormally bent).</p>
----------	---

Manual	4	<p>Cleaning</p> <p>Insert the products into an ultrasonic bath. Operate the ultrasonic device with tap water and 0.5 to 2.0% neodisher® MediZym for 10 to 30 minutes.</p> <p><i>Note: Follow instructions, observe water quality, concentrations and cleaning time specified by the manufacturer of the cleaning solution.</i></p>
---------------	----------	---

Manual	5	<p>Rinsing</p> <p>Rinse the products with tap water at 20-40°C for 1 min.</p> <p><i>Note: It is recommended to use deionized water.</i></p>
---------------	----------	--

Manual	6	<p>Drying</p> <p>Dry the products with compressed air until products are visibly dry.</p>
---------------	----------	--

Automatic	4	<p>Cleaning/Rinsing/Drying</p> <p>Place the instruments in the tray of the washer/disinfector's sliding trolley. Perform cleaning cycle with 0.2-1.0% neodisher® MediZym. Perform drying.</p> <p><i>Note:</i></p> <ul style="list-style-type: none"> Disinfection (thermal or chemical-thermal) is not required since the products are sterilized after cleaning. Follow instructions and concentrations specified by the manufacturer of the detergent solution. Follow the instructions of the washer-disinfector and verify that the success criteria have been met after each cycle as specified by the manufacturer. The final rinse step should be with deionized water. For other steps, follow the water quality defined by the manufacturer. Use only approved washers-disinfectors that comply with EN ISO 1588 in effect and are maintained and validated regularly.
------------------	----------	---

7	<p>Visual inspection</p> <p>Inspect the used products. Repeat steps 4-5-6 if the product is not visually clean or discard the product if it is damaged.</p>
----------	--

8	<p>Packaging</p> <p>Place the instruments in a paper-plastic pouch for steam sterilization in compliance in effect with ISO 11607 and EN 868 standards.</p> <p><i>Note:</i></p> <ul style="list-style-type: none"> For sharp devices that are not contained within a box, silicon tubes should be placed around the devices to prevent the packaging from being pierced. Seal the pouches as recommended by of the pouch manufacturer. If a thermo-sealer is used, the process must be validated, and the thermo-sealer must be calibrated and qualified.
----------	--

9	<p>Sterilization</p> <p>In accordance with the mandatory French regulations applicable in France (Prion Cycle):</p> <ul style="list-style-type: none"> Apparatus: class B Temperature: 134°C Time: 18 min Absolute pressure: 2.2 bar Minimum drying time: 20 min <p>In accordance with the current international standard ISO 17665-1, in effect:</p> <ul style="list-style-type: none"> Device: class B Minimum temperature: 132°C Minimum time: 3 min Absolute pressure: 2.2 bar Minimum drying time: 20 min <p>Check the physical-chemical indicators and cycle parameters.</p> <p><i>Note:</i></p> <ul style="list-style-type: none"> When sterilising multiple instruments in one autoclave cycle, ensure that the steriliser's maximum load is not exceeded. Place the pouches in the steam steriliser according to the recommendations given by the steriliser manufacturer. Use only a pre-vacuum air removal steam steriliser that meets the requirements in effect of EN 13060 (class B, small steriliser) and EN 285 (full-size steriliser), with saturated steam.
10	<p>Storage after sterilization</p> <p>Store the products in a dry, clean and dust-free environment at the temperature specified by the paper-plastic pouch by the steam sterilizer manufacturer.</p> <p><i>Note: Check the packaging and the medical devices before using them (packaging integrity, no humidity and use-by date). In case of damage, a complete rerun should be performed.</i></p>

11. Storage and transport conditions

Keep products away from direct sunlight and moisture.

12. Disposal

After use, in order to avoid any risk of infection and microbial risk, instruments must be placed in a secure container, used to collect cutting or sticking instruments (as needles or disposable bistouries) in accordance with good dentistry practices.

13. Symbols

	NiTi material		Quantity
	Stainless steel material		Assortment
	Manche «CLASSICS»		Catalogue number
	Non-sterile		Manufacturer
	Date of manufacture		Batch code number
	Can be sterilised in a steam steriliser (autoclave) at the specified temperature		

Year of EC marking: 2011
Version 4 – validation date: 03/2024

CFPM
2 bis, chemin du Loup 93290 Tremblay-en-France - FRANCE
Tel: +33 (0)1 45 91 31 82

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